


pyUSID & pycroscopy

Suhas Somnath, Chris R. Smith,
Nouamane Lanait, Stephen Jesse



ORNL is managed by UT-Battelle
for the US Department of Energy

Multitude of Instruments



Micro Raman Microscope



Atomic Force
Microscope (AFM)



AFM with Infrared
spectroscopy (AFM-IR)



Scanning
Tunneling
Microscope (STM)




Scanning
Transmission
Electron
Microscope (STEM)




AFM with Raman
spectroscopy

What we wanted



Instrument Tier




?



Interactive visualization, analysis on
supercomputers

Growing Data Sizes and Dimensionality

Evolution of Scanning Probe Microscopy Data



- Data sizes have grown from ~ 10 MB to ~ 1 TB in 10 years!
- Dimensionality ranges from 1D spectra to 7D hyperspectral datasets
- Cannot use laptops to analyze data


Instrumentation Software Inadequate for Analysis



- Software provided for controlling instruments typically only comes with basic data analysis capabilities
- Integrating user-developed functionality often impossible



Multitude of File Formats



Disjoint & Unorganized Communities



- Clustering
- Fit spectra ...



- Filter Image
- Register Image ...



- Fit Spectra
- SVD Filtering ...

- FFT Filtering
- SVD Filtering ...



- FFT Filtering
- Classify Images ...



- Register Images
- Clustering



MATLAB

Cannot Share Code Efficiently

- HIGHLY instrument-specific code
- Different programming languages
- Often licensed / costly software like Matlab
- Most popular sharing method = email!
- No centralized repository

Problems Opportunities in Imaging

- 1.** Closed science
 - a.** No traceability for data analysis
 - b.** Results not (readily) reproducible
- 2.** Multiple, incompatible, proprietary data formats
- 3.** Disorganized and unorganized communities
- 4.** No proper analysis software
- 5.** Growing data volumes, variety, and dimensionality

The Solution




Instrument Tier



Automated, standardized, modularized data acquisition



Instrument-agnostic, self-describing, model in HPC-friendly file format




Centralized, open repository for data processing

Expectation of Data Model

- Accommodate data of different shapes, dimensionalities, precision and sizes.
- Accommodate data without N-dimensional form
 - Compressed sensing / sparse sampling
 - Not all combinations of spectroscopic variables
 - Incomplete experimental data


Universal Spectroscopic & Imaging Data (USID)

- Data stored as 2D matrix of (position x spectral values) regardless of dimensionality
 - Ancillary datasets explain the data




USID – 1D spectra

Original N-dimensional form




Shape: (7,)
Quantity: Amplitude
Units: V

USID 2-dimensional form



USID – 2D Image



Original
N-D
form

Quantity: Intensity
Units: arb. units

=


Spectroscopic Position	P. I.	P. V.
0	0 0	-250 0
1	1 0	-125 0
2	2 0	0 0
3	3 0	125 0
4	4 0	250 0
0	0 1	-250 3.5
1	1 1	-125 3.5
2	2 1	0 3.5
3	3 1	125 3.5
4	4 1	250 3.5
0	0 2	-250 7
1	1 2	-125 7
2	2 2	0 7
3	3 2	125 7
4	4 2	250 7

S. I. arb.
S. V. arb.

Shape: (15, 1)

USID – Spectra on Grid (3D)


Original N-dimensional form



Quantity: Amplitude
Units: V

Suhas Somnath, somnaths@ornl.gov



USID 2-dimensional Form



S. I. Frequency



S. V. Frequency

USID – Images on a grid (4D)



?


USID – Images on a grid (4D)



Flatten
positions



Flatten
images



S.
V.

S. I.

-250	0	3.5
-125	0	3.5
0	0	3.5
125	0	3.5
250	0	3.5
-250	7	7
-125	7	7
0	7	7
125	7	7
250	7	7

P. I.


0	0
1	0
2	0
3	0
4	0
0	1
1	1
2	1
3	1
4	1
0	2
1	2
2	2
3	2
4	2

P. V.

-250	0
-125	0
0	0
125	0
250	0
-250	3.5
-125	3.5
0	3.5
125	3.5
250	3.5
-250	7
-125	7
0	7
125	7
250	7

USID - Instrument Agnostic Code

- Instrument-agnostic data allows instrument-agnostic code
- Single version of analysis and processing routine
- Brings multiple scientific communities together





Expectation of File Format

- Established standard in scientific research
- Store multiple datasets of different shapes, dimensionalities, precision and sizes.
- Scale very efficiently from few kilobytes to several terabytes
- Able to read and write data using any programming language including Python, R, Matlab, C/C++, Java, Fortran, Igor Pro, etc.
 - (without requiring installation of modules that are hard to install)
- Store metadata - experimental or analysis parameters
- Highly flexible and poses minimal restrictions on how the data can and should be stored.
- Compatible with cloud and high-performance computing paradigms (support parallel read and write)


Hierarchical Data Format (HDF5)

- A HDF5 file is a smart container
 - Capable of storing multidimensional datasets, Images, text, measurement parameters, etc.
 - Contents organized like traditional folders and files
 - **Groups** - Analogous to file folders
 - **Dataset** – 1 to N dimensional data
 - Integer, floating point, complex numbers etc
 - **Attributes** – {Key : value} pairs useful for describing data and experimental parameters, etc.
- Easily accessible – C, C++, python, Java....
- Tree structure + nomenclature +attributes are **records of workflow** applied to dataset
- Parallel read / write, HPC & cloud compatible



Traceability and Reproducibility


Raw, intermediate, final results stored in same file, leveraging hierarchy



All measurement and analyses parameters stored for repeatability

```
Measurement_000 (2088, 4)
Group size = 1
Number of attributes = 43
BE_actual_scan_time_[s] = 0.004
BE_amplitude_[V] = 4
BE_auto_smoothing = auto smoothing on
BE_band_edge_trim = 0.14614
BE_band_smoothing_[Hz] = 2560.5
BE_band_width_[Hz] = 60000
BE_bins_per_band = 0
BE_center_frequency_[Hz] = 370000
BE_phase_variation = 1
BE_points_per_BE_wave = 0
BE_repeats = 8
BE_signal_type = chirp-sinc hybrid
BE_time/pixel_[s] = 0.004
File_date_and_time = 06-Feb-2015 11:46:13
File_file_name = BELine
File_file_path = C:\Users\Administrator\Documents\Asylum Research Data\150205\
File_file_suffix = 9
IO_AO_amplifier = 1
IO_AO_range_[V] = +/- 10
IO_Analog_Input_1 = +/- 1V, FFT
IO_Analog_Input_2 = +/- 10V, mean
IO_Analog_Input_3 = off
IO_Analog_Input_4 = off
IO_DAQ_platform = NI 5412/5122
IO_rate_[Hz] = 4000000
data_type = BELineData
grid_contact_set_point_[V] = 1
grid_current_row = 1
grid_cycle_time_[s] = 0.05
grid_lift_height_[m] = 5.0E-8
grid_nap_mode = nap mode off
grid_num_cols = 128
grid_num_rows = 128
grid_scan_time_/_line_[s] = 1
grid_time_remaining_[h:m:s] = 10
grid_total_time_[h:m:s] = 10
machine_id = mac109728.ornl.gov
num_bins = 44
num_pix = 16384
num_udvs_steps = 1
platform = Darwin-17.4.0-x86_64-i386-64bit
pycroscopy_version = 0.60.orc1
timestamp = 2018_05_03-15_23_37
```

Facilitating Data Mining



ORNL Data Catalog

Datasets Organizations Groups About

Home / Datasets

Search datasets...

Order by: Relevance

10 datasets found

Formats: h5 x

BFO grain boundary 4D STEM
Highly strained polymorph bismuth ferrite (BFO) thin film.
h5

BFO S-PFM
Single frequency piezoresponse force microscopy scan of Bismuth ferrite.
h5

BEPS on PMN-PT
Band Excitation Polarization Switching on (100)-textured 0.70PbMg2/3Nb1/3O3-0.3PbTiO3 thin films.
h5

STS Data La58Ca38MnO3
Atomic-scale surface structure of epitaxial La58Ca38MnO3 thin film obtained from scanning tunneling spectroscopy.
h5

BELLine Imaging of BFO
Band Excitation Imaging of Bismuth ferrite obtained from an Asylum Research AFM

Datasets Organizations Groups About

Organizations

Materials Science a... 10

Groups

There are no Groups that match this search

Tags

AFM 3

BFO 2

imaging 2

microscopy 2

spectroscopy 2

SPM 2

visualization 2

atmoic-force-micros... 1

BEPS 1

epitaxial-strain-re... 1

Show More Tags



Formats

h5 10

Licenses

(Purely) Programmer-Driven Solution

Software connecting scientific instruments to supercomputers



- **Successes:**
 - Easy to use – Point-click
 - Fast – on super-computer
- **Shortcomings:**
 - Very long development cycle
 - Very expensive
 - Brittle (points of failure)
 - Scientists had no control!!

Expectation from Software

- Easy to learn and understand
- Strong support-base
- Established community standard
- Straightforward to implement and maintain
- Optimized libraries for scientific and numeric algorithms
- Access to existing imaging related packages
- Affordable
- Scalable to multiple CPU cores + distributed computing


Two-prong (inverse) Software Strategy

Mode	Exploratory / Development	Production
Compute Resource	Instrument-local Workstation 	Cloud compute 
Who?	Domain scientists	Collaboration with computational facility
Algorithms	Simpler but approximate	Complex but accurate


Python for Scientific Research

Very easy to learn + code

Numerous, **powerful** libraries for science



NumPy




- Facilitates innovation
 - More robust code
 - Improved adoption of new methods / standards
- Accelerates scientific progress

Cross-platform



scalable



Established standard for:

- | | |
|---|---|
| <ul style="list-style-type: none">• Microscopy• Microbiology• Deep learning | <ul style="list-style-type: none">• Data science• Neutron science• More...! |
|---|---|

Strong user community



All for a princely sum of **\$0!**

See Jake Vanderplas' Pycon 2017 Keynote talk:


<http://www.youtube.com/watch?v=ZyjCqQEUA8o&t=19m20s>

Two Software Packages

- Written in Python
- Open source & free
- Written by scientists
- Data centric
- Instrument-independent data model in HDF5
- Instrument-independent analysis algorithms
 - Reusable across scientific domains



2,800+ downloads



121,000+ downloads

pyUSID – Software Organization

pyUSID

io

- Tools to read and write to USID HDF5 files
- *USIDataset* – slicing, visualization, reduction, reshaping of 2D USIDsets
- Data Translators
 - Images & numpy arrays to HDF5


Processing

- *Process* –
 - Frame science as compute problem
 - Parallel & scalable computation
 - Multi-core (joblib)
 - Multi-node (MPI)
 - Piecewise, resume
 - Manages memory

Visualization

- Plotting utilities
- Interactive jupyter widgets for N-dim datasets

Scaling on HPC



* Preliminary results on ORNL's CADES SHPC Condo

** Lower than ideal trend because compute job was too small to use many nodes

Pycroscopy - Software Organization

pycroscopy

io

- Data translators (proprietary formats to HDF5)

Visualization

- Science-specific plotting & interactive jupyter widgets

Simulation

- AFM Force-distance ...

Analysis

- Physical model specific
- Fitting to model, etc.

Processing

- Physical model agnostic
- Image filtering, registration,
- Multivariate analysis

Software Organization

Science / data
analytics
applications


File &
computing
tools




Well documented

Beginner topics


To learn how to use pyUSID, Please go through the following documents in the recommended order:




01. Primer to HDF5 and h5py




02. The USIDataset



03. Translation and the NumpyTranslator



04. Plotting utilities



05. Utilities for reading h5USID files

To view the file ('ifft2'). Remember it is necessary to use the inverse transform. Also the inverse transform is symmetric about the center, resulting in the inverse being times smaller than the original kept.

image_filter
image_filter

```
fig, axes =  
    for axis, img, title in zip(axes, [image_raw, image_filtered], ['original', 'filtered']):  
        _ = px.plot_utils.plot_map(axis, img, cmap=plt.cm.inferno,  
                                  x_size=x_edge_length, y_size=y_edge_length, num_ticks=5)  
        axis.set_title(title)  
fig.tight_layout()
```

`reshape_to_Ndims(h5_main, h5_pos=None, h5_spec=None, get_labels=False, verbose=False, sort_dims=False)` [\[source\]](#)

Reshape the input 2D matrix to be N-dimensions based on the position and spectroscopic datasets.

- Parameters:
- **h5_main (HDF5 Dataset)** – 2D data to be reshaped
 - **h5_pos (HDF5 Dataset, optional)** – Position indices corresponding to rows in *h5_main*
 - **h5_spec (HDF5 Dataset, optional)** – Spectroscopic indices corresponding to columns in *h5_main*
 - **get_labels (bool, optional)** – Whether or not to return the dimension labels. Default False
 - **verbose (bool, optional)** – Whether or not to print debugging statements
 - **sort_dims (bool)** – If True, the data is sorted so that the dimensions are in order from fastest to slowest. If False, the data is kept in the original order. If *get_labels* is also True, the labels are sorted as well.

Returns:


- **ds_Nd (N-D numpy array)** – N dimensional numpy array arranged as [positions slowest to fastest, spectroscopic slowest to fastest]

Intermediate topics


To learn how to write to h5USID files, write data processing classes, or adding functionality to pyUSID, go through these additional documents in the recommended order: Those interested in contributing to pyUSID are encouraged to read our [guidelines for contributing code](#)




06. Utilities for handling



07. Speed up




08. Utilities that assist



Entering the USID Ecosystem

- hdf5 file is the hub for all operations
- Analysis, processing, visualization available after translation to .hdf5



Jupyter Notebooks



Jupyter Notebook

A screenshot of a Jupyter Notebook interface. On the left, there's a sidebar with the 'jupyter' logo, a 'Welcome to the' message, and a 'Run some Python' section. The main area shows a notebook cell titled 'Exploring the Lorenz System'. It contains text about the Lorenz system, differential equations, and atmospheric convection. Below that is another cell with code for interactively exploring the system using sliders for parameters α , β , and ρ . The final part of the screenshot shows a 3D plot of the Lorenz attractor, which is a complex, chaotic attractor with two symmetric butterfly-shaped lobes.

- Interactive documents
- Exploratory programming
- Code
- Text
- Images
- Interactive – slice through data, pan, move, rotate ...

Towards Open & Reproducible Science

Aim – All scientific journal papers accompanied with:

- Jupyter notebook that shows all analysis (raw data → figures).
- Data with DOI number
- Containers – software environment

The screenshot shows a Nature Communications article page. At the top left is the journal logo. Below it, the word "ARTICLE" is in blue. To its right is an orange "OPEN" button. Underneath is the DOI: 10.1038/s41467-017-02455-7. The title of the article is "Ultrafast current imaging by Bayesian inversion". Below the title is a list of authors: S. Somnath^{1,2}, K.J.H. Law^{1,3}, A.N. Morozovska^{1,4}, P. Maksymovych^{1,2}, Y. Kim⁵, X. Lu^{1,6}, M. Alexe⁷, R. Archibald^{1,3}, S.V. Kalinin^{1,2}, S. Jesse^{1,2} & R.K. Vasudevan^{1,2}. The abstract begins with: "Spectroscopic measurements of current-voltage curves in scanning probe microscopy is the earliest and one of the most common methods for characterizing local energy-dependent electronic properties, providing insight into superconductive, semiconductor, and memristive behaviors. However, the quasistatic nature of these measurements renders them extremely slow. Here, we demonstrate a fundamentally new approach for dynamic spectroscopic current imaging via full information capture and Bayesian inference. This general-mode I-V method allows three orders of magnitude faster measurement rates than presently possible. The technique is demonstrated by acquiring I-V curves in ferroelectric nanocapacitors, yielding >100,000 I-V curves in <20 min. This allows detection of switching currents in the nanoscale capacitors, as well as determination of the dielectric constant. These experiments show the potential for the use of full information capture and Bayesian inference toward extracting physics from rapid I-V measurements, and can be used for transport measurements in both atomic force and scanning tunneling microscopy."

Suhas Somnath, somnaths@ornl.gov


Jupyter notebook associated with paper

Visualizing Filtering Results

Run the next cell to see what the IV curves look like for the chosen filter parameters. If the current set of filter parameters do not work as well, change the parameters in the previous cell and try again.

Once the results look good, proceed to filter the entire dataset

```
In [7]: # Test filter on a single line:  
row_ind = 50  
filt_line, fig_filt, axes_filt = px.processing.gmode_utils.test_filter(h5_main[row_ind], filte  
r_parms, samp_rate,  
                                         show_plots=True, use_rai  
nbrow_plots=False)  
fig_filt.savefig(os.path.join(other_figures_folder, 'FFT_filter_on_line_{}.png'.format(row_in  
)), format='png', dpi=300)  
  
raw_row = np.reshape(h5_main[row_ind], (-1, pts_per_cycle))  
filt_row = filt_line.reshape(-1, pts_per_cycle)  
  
fig, axes = px.plot_utils.plot_loops(single_A0, [raw_row, filt_row], dataset_names=['Raw', 'Fi  
ltered'],  
                                         line_colors=['r', 'b'], x_label='Bias (V)', title='FFT Fi  
ltering',  
                                         plots_on_side=3, y_label='Current (nA)',  
                                         subtitles='Row: ' + str(row_ind) + ' Col: ')  
fig.savefig(os.path.join(other_figures_folder, 'Example_filtered_loops_from_line_{}.png'.forma  
t(row_ind)), format='png', dpi=300)
```



DOI associated with data (raw → paper figures)

Oak Ridge National Laboratory 10.13139/OLCF/1410993 [Download](#)

Authors

Somnath, Suhas	somnaths@ornl.gov
Law, Kody	lawkj@ornl.gov
Morozovska, Anna	anna.n.morozovska@gmail.com
Maksymovych, Petro	maksymovychp@ornl.gov
Kim, Yunseok	ykim943@gmail.com
Lu, Xiaoli	xliu@xidian.edu.cn
Alexe, Marin	M.Alexe@warwick.ac.uk

Pycroscopy - Supporting User Research

Before 2016


```
clear all
%
% Mycroscopy data is 4D in which each (x,y) position in a 2D scan
% contains a 2D image of scanned electron intensity. The data is read
% in as a series of 2D images of binning factors (2x2x2x2). The data has to be
% hand sorted into a 4D array.

% Define file path
file_path = 'C:\Users\1\desktop\Mycroscopy\slicing_and_dicing'


% determine all the file names of image files within the folder
% d = dir(file_path); use the 'd' command to retrieve all file names. Then
% for all files in the folder, read them in as 2D images
N_ronchigrams = n-2; Masterlist the number of images in the folder

% some additional information that will be needed when loading the data
size_ronc = 256; Number of pixels along one edge of ronchigram
size_scan = 128; Number of pixels along one edge of scanned image. Number of pixels along one edge of ronchigrams
size_ronc/binning_factor = 128; % because the data is so large, it will be beneficial
% perform batch loading of all of the files in the folder
ronc_mean = zeros(N_ronchigrams, size_ronc/binning_factor, size_ronc/binning_factor);
plot_cond = 1; Turn plotting on (1) or off (0)

for k1 = 1 : N_ronchigrams
    fraction_loaded = k1/N_ronchigrams;
    disp(fraction_loaded)
    file_name = d(k1).name;
    get_file_name_for_ronc_mean = imread(fullfile(file_path, file_name));
    ronc_mean(k1) = double(ronc_mean);
    % convert type
    if(plot_cond==1)
        figure(1)
    end
end
```



Since 2016



Scripts + complicated, Matlab GUI

Set of simple Jupyter notebooks

Written by dedicated software engineer

Written by material scientists

Not customizable

Completely customizable.

2-3 hours of training before use

Notebooks include instructions. NO training required!

Deployed only on two offline workstations due to licensing restrictions = queue

Each user gets VMs with jupyter notebook server

Will remain on off-line desktops

In the process of switching to computations on clusters

Outreach - workshops

MM_2018_Workshop

Tutorial materials for the Microscopy and Microanalysis Workshop

Prerequisites:

Please follow [instructions here](#) to install Anaconda and pycroscopy

More information on the contents:

Please see our website. You should be able to find all information ranging from basic concepts to advanced topics such as machine learning and deep learning enabled by pycroscopy, reading different microscopy file formats, tutorials on

Getting in touch

Time	Details
8:00-9:00 AM	Setup and installation
9:00-9:30 AM	Introduction to machine learning in python
9:30-10:15 AM	Linear unmixing methods for spectral analysis
10:15-10:30AM	Coffee Break
10:30-12:00 PM	Nonlinear unmixing: overview and applications
12:00-1:00 PM	Lunch on your own
1:00-2:00 PM	Traditional classification methods
2:00-2:45 PM	Clustering methods applied to 4D STEM data
2:45-3:00 PM	Break
3:00-4:30 PM	Neural networks and deep learning for microscopy images
4:30-5:00 PM	Wrap up/Discussion/Summary

Imaging and Spectral Data Analysis in Python

We present a set of jupyter notebooks and presentations that go over the basics of:

- Image and spectral data analysis in Python
- [USID](#), [pyUSID](#), and [pycroscopy](#) as tools for analyzing large and multidimensional imaging / spectroscopy datasets

This material was developed as part of the [2018 CNMS User Meeting](#)

See our [list of tutorials](#) for more information


Getting started

(Recomm) [CNMS_ML_in_MS_Workshop_2018](#)


1. (one tip)
 - i. Fc
 - ii. AlCNMS User Meeting Workshop 2018: Machine learning in Materials Science with Python
Welcome! This workshop will go over machine learning methods for materials science with python. The full agenda is published below. For remote participants, please click <https://bluejeans.com/1766282540>
2. (one tip)
 - i. cli
 - ii. Cl
 - iii. Do
 - iv. UrThe workshop will be conducted on August 15, 2018 at ORNL. The livestream begins at 9AM EDT.
To get the most out of this workshop, you may want to work along with us on your own laptop. If you do, we highly recommend that you:
 1. Download and install the following:
 - Anaconda5.2 for python 3.6 - <https://www.anaconda.com/download/>
 - Pycroscopy- <https://pycroscopy.github.io/pycroscopy/install.html>
 - HDFview- <https://support.hdfgroup.org/products/java/release/download.html>
 2. Familiarize yourself with basics concepts.Links to quick and helpful tutorials are available here - https://pycroscopy.github.io/pyUSID/external_guides.html.
Topics include:
 - Basic python usage
 - Jupyter notebooks for running code
 - Numpy for numeric operations

Pycroscopy - Scientific Advancements


3,500x faster imaging via adaptive signal filtering, linear unmixing of signals




200x faster spectroscopy via Bayesian inference




Separating uncorrelated data from correlated data to clean images



Identifying invisible patterns using multivariate analysis




Simplified navigation multidimensional data - users




Software Progression

Scaling up Computing:



New Scientific Domains:



Lots more information available online

<https://pycroscopy.github.io/pyUSID/about.html>

The screenshot shows the pyUSID documentation website. The top bar includes a logo, the version '0.0.5', and a search bar. The sidebar on the left lists navigation links such as 'What?', 'Why?', 'Who?', 'Getting Started', 'Installation', 'Tutorials on Basics', 'Package Organization', 'Examples & Tutorials', 'Data Schema and File Format', 'Guidelines for Contribution', 'Frequently asked questions', 'Contact us', 'Credits', 'What's New', and 'Upgrading from Matlab'. The main content area displays the 'pyUSID' page, which includes a brief description of the Python framework for storing and processing spectroscopy, imaging, observational / experimental data.

pyUSID

Python framework for storing and processing spectroscopy, imaging, observational / experimental data

What?

- The [Universal Spectroscopy and Imaging Data \(USID\) model](#):
 - facilitates the representation of spectroscopic or imaging data regardless of its origin, modality, size, or dimensionality.
 - enables the development of instrument- and modality- agnostic data processing and analysis algorithms.
 - is just a definition or a standard, something tangible and concrete.
- pyUSID is a [python](#) package that provides three pieces of functionality:
 1. **io**: Primarily, it enables the reading and writing of USID in hierarchical JSON files (referred to as h5 files).
 2. **viz**: It has handy tools for visualizing general scientific data.

The screenshot shows the pycroscopy documentation website. The top bar includes a logo, the version '0.60.2', and a search bar. The sidebar on the left lists navigation links such as 'What?', 'Why?', 'How?', 'Who?', 'Papers / Conferences', 'Getting Started', 'Installation', 'Tutorials on Basics', and 'Package Organization'. The main content area displays the 'pycroscopy' page, which features a cartoon snake holding a graduation cap and a diploma, with the text 'pycroscopy' overlaid.

<https://pycroscopy.github.io/pycroscopy/about.html>

Pycroscopy

Scientific analysis of nanoscale materials imaging data

Note

We are hiring! We are looking for full-time python software developers who can support the scientists' needs at the Center for Nanophase Materials Science (CNMS) in addition to taking pycroscopy forward. More details available [here](#)

What?

- pycroscopy is a [python](#) package for processing, analyzing, and visualizing multidimensional imaging and spectroscopy data.
- pycroscopy uses the [Universal Spectroscopy and Imaging Data \(USID\) model](#) as its foundation, which:
 - facilitates the representation of any spectroscopic or imaging data regardless of its origin, modality, size, or dimensionality.
 - enables the development of instrument- and modality- agnostic data processing and analysis algorithms.
- pycroscopy uses a data-centric model wherein

Thank you

Questions?



Science, data analysis, reading proprietary instrument data -

<https://groups.google.com/forum/#!forum/pycroscopy>

File / data tools related- <https://groups.google.com/forum/#!forum/pyusid>

This presentation -

<https://docs.google.com/presentation/d/1K9SCTQU6KSSAewlkiwWz77FDfKSzZJqZPnTy2z9OujU/edit?usp=sharing>